



Labvill

Laboratory Milling to Predict Wheat Behavior at Plant Level

Arnaud Dubat

Business development Director



25/10/2016

Presentation Map





Presentation Map







Milling: a science made of several steps!



Wheat kernel

- ✓ Wheat kernel is not a rice kernel : because of the crease!
 - To remove the bran and obtain white flour, the miller needs to act progressively.
 - First step = breaking ("opening" the grain)
 - Grooved/corrugated rolls
 - Second step = reduction (avoiding breaking bran particles)
 - Smooth rolls

There are no mills in the world that are not using smooth rolls (common soft/hard wheat).





Wheat Milling economics



Extraction rate



How much could you save?

✓ Total yield : 98%

✓ Average flour yield : 78%

1% flour extraction extra means:

+1476 €/day

+500 000 €/year

and the set of the set

What is the goal of a laboratory mill?

In the milling industry, wheat is characterized by its milling behavior (resistance to crushing and extraction rate) and by the quality of the flour produced.



ZabMill is designed to anticipate, in the laboratory, the behavior of wheat in the industrial mill, and to evaluate its extraction potential and the quality of the flour produced.









*LabMill was developed within the Milling Performance Consortium (AFSA, Arvalis-instituted du végétal, ANMF, Danone Vitapole, INRA, IRTAC, Ulice, CHOPIN Technologies).



KPM





*LabMill was developed within the Milling perfromance Consortium (AFSA, Arvalis-institut du végétal, ANMF, Danone Vitapole, INRA, IRTAC, Ulice, CHOPIN Technologies).









*LabMill was developed within the Milling performance Consortium (AFSA, Arvalis-institut du végétal, ANMF, Danone Vitapole, INRA, IRTAC, Ulice, CHOPIN Technologies).









*LabMill was developed within the Milling performance Consortium (AFSA, Arvalis-institut du végétal, ANMF, Danone Vitapole, INRA, IRTAC, Ulice, CHOPIN Technologies).









*LabMill was developed within the Milling performance Consortium (AFSA, Arvalis-institut du végétal, ANMF, Danone Vitapole, INRA, IRTAC, Ulice, CHOPIN Technologies).









*LabMill was developed within the Milling performance Consortium (AFSA, Arvalis-institut du végétal, ANMF, Danone Vitapole, INRA, IRTAC, Ulice, CHOPIN Technologies).



 KPM
 Patent FR N° 0905572 | US N°9 067 210

CHOPIN

LabMill - Key features

Modern and user-friendly interface (7" touch screen)

✓ Very precise and adjustable feed system

- Feeding hoppers equipped with a precision scale
- Provides a **constant flow of grain**, adapted to every type of wheat

✓ Individually adjustable cylinders

Adjust settings to mimic industrial milling

✓ Centrifuge sifting

- Compact, economic and easy maintenance
- Easy to change sieves
- Time of sifting adjustable

Non automatic mill

• Because millers want to see the different products









LabMill provides high extraction rates...



- Range: 66.2%
 to 81.5%
- Average hard wheat: 77%
- Average soft wheat: 75%

 Note: samples tempered to 16.5% moisture content for 24 hours





LabMill shows the great differences on wheat behavior during milling

- **150 wheat samples**, coming from **14 countries** distributed on the **5 continents**
- LabMill with the same protocol -> change behavior and results if we adjust the LabMill
- Flour production is highly variable from one sample to other.
- Not only the **global quantity of flour varies** (extraction rate changes from 59.4% to 77.5%)







LabMill shows the great differences on wheat behavior during milling

But also the distribution among the different milling steps :

- breaking step produces between 14%
 and 58% of the total flour
- **sizing** step between 13% and 23%
- **reduction** step between 27% and 70%





...and low ash and damaged starch contents!



Range: 0.42% to 0.76%
Average: 0.58%
This corresponds to :
French : Type 55
American : All purpose & bread
Argentinian : "000"
German : 550

Variety	SOISSONS	ΡΑΚΙΤΟ	APACHE	ARKEOS1	ARKEOS2	CAPHORN	CALCIO	CROUSTY	MATHEO
Ash (%)	0.57	0.55	0.56	0.45	0.5	0.62	0.59	0.51	0.66
Damaged starch (UCD)	16.63	15.38	13.53	6.77	10.03	20.28	15.52	7.75	14.8

Milling Quality Consortium – November 4, 2015

Presentation Map







Experimental setup

- 150 wheats from soft to hard originating from 17 countries.
- 2 LabMill (A & B)
- 3 analysts
- Observation period 4 month (April-August 2016)
- Each sample (1700g) was tested 5 times on each mill.
 - Repeatability is measured as the average standard deviation of the 5 subsamples tested on the 2 mills $(Sr_{LabmillA} + Sr_{LabmillB})/2$
 - Reproducibility is measured as the standard deviation of the 10 subsamples tested on the 2 mills (A and B)





Example:1st Break Flour



Our observations cover a wide range of B1 flour [3%-25%) In this range repeatability and reproducibility remain constants

	Average CVr(%)**	Average CVR(%)**
Extraction rate	1,9%	4,2%

October 8, 2018





Summary

Parameter	Range (%)	Average	Repeatability CV% Calculated on the mass of product (in g)	Reproduciblity CV% Calculated on the mass of product (in g)
Extraction rate	56-76	70,3%	1,2%	1,7%
1st Break Flour	3-25	10,1%	1,9%	4,2%
2 nd Break Flour	6-16	11,2%	1,4%	4,0%
Sizing Flour	8,5-18,5	12,8%	2,3%	4,2%
1st Reduction flour	8-35	24,5%	4,7%	6,6%
2 nd Reduction flour	5-14	7,8%	7,8%	10,7%
3rd Reduction flour	3-7	3,2%	10,8%	13,4%
Coarse Bran	8-24	14,4%	3,3%	6,1%
Fine Bran	4-12	6,5%	2,7%	8,2%
Shorts	3-18	7,8%	6,3%	12,2%





Conclusions

- This first approach shows
 - A very good repeatability, higher CVs are related to streams with relatively little amount of product, then a minimum difference in weight has higher impact on CVs.
 - Reproducibility is also good taking into account the large variability of test conditions.
 - Last but not least the LabMill shows to perform exactly the same on all kind of wheat making it a mill suitable whatever the wheat hardness is.





Presentation Map







Origin of the idea

 150 wheat, submitted to the same milling, have shown extremely diverse behaviors.

• Based on these observations, could we develop a way to easily characterize wheat behavior during milling?

• What could be the benefit of such measurement?





Resistance and dissociation Indexes

 Based on these observations we can characterize wheat behavior during milling with a simple Resistance to Crushing and Dissociation Indexes





The Reduction index







Milling Performance Index (MPi), describes the wheat behavior during milling.



October 8, 2018

Practical example 1 5 samples with same Extraction rate and different MPi

MPi	Extraction rate (%)	Meal/B1 (%)	Large Middlings B1 (%)	Fine Middlings B1 (%)	Flour B1 (%)	Coarse Bran (%)	Large Middlings/B2 (%)	Fine Middlings/B2 (%)	Flour B2 (%)	Fine Bran (%)	Fine Middlings/Sizing (%)	Flour Sizing (%)	Fine Middlings/Red1 (%)	Flour/Red1 (%)	Fine Middlings/Red2 (%)	Flour/Red2 (%)	Fine Middlings/Red3 (%)	Flour/Red3 (%)	Feed material → Sizing (%)	Feed Material → Reduction (%)
233	70,0	57,6	14,3	12,1	15,9	17,1	7,9	17,9	14,5	5,2	5,3	11,4	16,7	18,1	10,1	6,7	7,4	2,8	22,1	35,3
333	70,0	58,6	16,0	11,8	13,5	15,4	9,1	21,3	12,6	5,7	6,8	12,4	19,0	20,5	12,0	6,9	8,7	3,4	25,0	39,8
444	69,8	61,5	15,4	10,9	12,0	16,5	10,1	22,6	12,0	5,6	7,9	11,9	18,9	22,0	11,2	7,7	7,9	3,4	25,5	41,3
555	70,4	64,5	18,1	9,5	7,9	13,1	13,9	26,8	10,4	7,2	11,8	12,8	23,6	24,2	13,4	10,2	9,1	4,4	32,0	48,2
644	70,4	66,0	16,9	9,6	7,4	12,7	13,9	28,2	11,0	7,4	10,2	13,1	21,7	25,9	13,3	8,4	9,6	3,7	30,8	48

All % are given as a % of entering wheat



Practical example 1 5 Samples with same Extraction rate and different MPi

Looking at	Average	Range	Variation
Feeding 2nd Break	6 t/h	1t/h	+/- 17 %
Feeding Sizing	2,5t/h	1t/h	+/- 40 %
Feeding Reduction	4,2 t/h	1,3 t/h	+/- 31 %

- Same flour extraction can ne obtained with samples responding very differently to the milling process.
- In the industry this can generate difficulties at different stages (grinding, sifting...)



-eed numbers in t/h are hypothetical, helping comparison --not example of existing situation

Practical example 2 5 samples with same MPi and different Extraction rate

	MPi	Extraction rate (%)	Meal/B1 (%)	Large Middlings B1 (%)	Fine Middlings B1 (%)	Flour B1 (%)	Coarse Bran (%)	Large Middlings/B2 (%)	Fine Middlings/B2 (%)	Flour B2 (%)	Fine Bran (%)	Fine Middlings/Sizing (%)	Flour Sizing (%)	Fine Middlings/Red1 (%)	Flour/Red1 (%)	Fine Middlings/Red2 (%)	Flour/Red2 (%)	Fine Middlings/Red3 (%)	Flour/Red3 (%)	Feed material → Sizing (%)	Feed Material $ ightarrow$ Reduction (%)
	344	70,4	59,9	16,1	11,4	12,4	15,7	9,8	22,2	11,9	5,8	7,7	12,3	19,5	21,3	11,5	8,0	8,2	3,4	25,9	41,3
	344	71,7	58,4	16,4	11,8	13,4	14,5	9,4	21,8	12,3	5,9	6,9	12,9	18,5	21,6	10,7	7,8	7,7	3,0	25,8	40,5
1	344	72,9	59,6	16,2	11,0	12,8	15,4	9,9	21,8	12,2	5,1	7,2	13,5	16,0	23,8	9,0	7,0	6,3	2,7	26,0	40,1
「いうて	344	73,8	57,3	17,4	12,0	13,2	14,2	8,9	21,5	12,5	5,8	6,8	13,6	15,7	24,2	8,5	7,2	5,9	2,6	26,3	40,2
N.	344	74,8	58,8	17,4	11,3	12,4	13,0	10,2	23,6	11,8	5,6	8,4	13,3	16,3	26,5	8,9	7,4	6,3	2,6	27,5	43,3
Mar Andrew	344 All % ar	74,8 e given a	58,8 as a % c	17,4 of ente	11,3 ring wh	12,4 eat	13,0	10,2	23,6	11,8	5,6	8,4	13,3	16,3	26,5	8,9	7,4	6,3	2,6	27,5	



Practical example 2 5 Samples with same MPi and different Extraction rate

Looking at	Average	Range	Variation
Feeding 2nd Break	5,9 t/h	0,2 t/h	+/- 3 %
Feeding Sizing	2,7 t/h	0,2 t/h	+/- 7 %
Feeding Reduction	4,15 t/h	0,3 t/h	+/- 7 %

- Because it is based on the observation of wheat behavior during crushing MPi classifies wheat according to their response to milling stress.
- This allows a new and useful way to look at Laboratory milling results.



Feed numbers in t/h are hypothetical, helping comparison -

not example of existing situation



Wheat milling behavior at a glance : and a graphical visualization of the results





Let's start REALLY analyzing wheat milling

• The Milling Value index is a new parameter that helps understanding wheat behavior during milling.

• Certain wheat can have high extraction potential but grinding characteristics not suitable for certain mil diagrams.

 Selecting wheat according to their potential extraction AND behavior during milling makes sense to ensure that incoming wheat is adapted to millers needs.





Presentation Map







Comparison with the French reference **Pilot mill** (ENILIA-ENSMIC, Surgères) ARVALIS Institut du végétal







CHOPIN Technologies

introduction

• The objective is to compare milling behavior of the CHOPIN LabMill with the French reference Pilot mill on 22 wheat samples.

• The LabMill was used according to the classical tempering protocol (ISO 27971) 16% moisture, 24H resting time. Only 2 reductions were used instead of 3. (Tess done by Arvalis-institut du végétal)

• At the Pilot mill level, all tempering were conducted at 16% moisture (except for 4 wheat tempered at 16,5%). Resting time was 48 hours. Tests done by Enilia-Ensmic.





What are we comparing??











Br 1



caisse



CL3			2*1120	germes
Détacheur			3*425	B4
Ecartement	1,67/100		5*112	
Vitesse L		F1		
Vitesse R			3*250	CL4/C4
			C2	





F1

,67/100

C1

Detachour

/itesse

cartiement

Offérentie

		1/2 c	aisse	
C3+filtre			3*132	CL4/C4
Détacheur		F1		
Ecartement	contact			
Vitesse				
Différentiel				

		1/2 caisse	
C5		1*112	
Détacheur			
Ecartement	contact	F1	
Vitesse		3*132	Rem bis
Différentiel	F	Rem blanc	



Diagram LabMill and products

Reduction 1

Flour

Reduction 2

Flour

KPM ANALYTICS October 8, 2018

0,75 mm son Pis C3

A direct prediction of B1 flour !

Labmill Predictor For industrial B1 flour	Coef Correlation (r)				
% Farine B1	0,964				
% Total middlings	-0,910				
% Fine middlings B1	0,891				
% Fine middling total	-0,889				
% Fine middling B2	-0,878				
% Flour C2	-0,861				
% Fine middling Sizing	-0,854				

Wheat behavior at B1, directly related to their resistance to crushing, confirms the strong prediction potential of the LabMill. This results is important as it shows the achievable analogy between Laboratory and industrial milling.







Prediction of main flour production zone

Sizing flour % Prediction

Break flour % Prediction

32,00%

34,00%

36,00%

38,00%

Prédiction %Farine Convertissag

40,00%



42,00%





Bran & shorts prediction

Shorts % Prediction

% remoulage total surgères = 0,000437 + 0,9991 Prédiction %Remoulages







BRAN % Prediction

% sons Surgères = 0,00001 + 1,000 Prédiction % Sons

Extraction rate prediction





E KPM Octo

Conclusion

• Results analysis shows a great potential of the LabMill to predict Surgères Pilot mills results.

 At this stage, and taking into consideration that 1/Laboratory and industrial diagram are strongly different and 2/normal measurement uncertainty at lab and on line, these results are very encouraging and show a great potential to create the missing link between the Laboratory and the production plant.





Presentation Map





How to measure blending law?

Analytics

- 2 wheats : Hard and Soft
- Blending :
 - 100% Hard
 - 75% Hard + 25% Soft
 - 50% Hard + 50% Soft
 - 25% Hard +75% Soft
 - 100% Soft
- Tempering 24h, 16% H2O
- LabMill (standard settings)
- Analytical
 - Alveograph HC
 - Spectralab (Moisture, protein, ash)
 - SDmatic (Starch Damage)
 - CHOPIN-SRC

Data study

- Calculate blending values according to the wheat blending (calculated values)
 - Example 50/50
 - Calculated value = 50% real Hard wheat value + 50% real Soft wheat value
- Compare calculated value vs real value
- Assess coefficient correlation (r²)





MPI based indexes



- It is possible to base wheat blending on
 - Resistance Index (the most important of the 3 indexes) (r²=0,94)
 - Apparent Hardness index (r²=0,99)





Break Flour



 1st and 2nd Break flour perfectly respond to blending law (r²=0,99 and 0,84 respectively)





Sizing flour



- Sizing flour does not look to respond blending law.
- This is directly related to very little variation in flour production at this stage between soft and hard wheat (0,5%)
- In this study, this value is a constant.

	Measured	Calculated
Average	13,3%	13,6%
Mini	12,8%	13,4%
Maxi	13,9%	13,9%



October 8, 2018

Conversion flour



- Red 1 responds well to blending law (r²=0,84)
- Red2 responds less (r²=0,60) but, here again, variation between hard and soft is very limited (2%) (same for Red3 –not shown here)





Alveograph







- Alveograph data respond very well to blending law. (R²= 0,98 for W and P, 0,89 for Ie)
- This was already known blending flour, it is confirmed blending wheat on the LabMill



Conclusions

- Using the blending law with the LabMill works* for predicting:
 - MPI and apparent hardness values
 - Different streams products
 - Flours quality
- This means that the LabMill is a suitable tool to:
 - Predict the final result and optimize the blending of 2 wheats
 - Calculate the « ideal » wheat to correct a deficient one.
 - Example : you have a wheat with apparent hardness 80, you want to blend it 50/50 with another one and have final apparent hardness 85
 - 85 = 0,5*80 + 0,5*? → ? = (85-(0,5*80))/0,5 = 90





Conclusions

- Laboratory milling is more than just grinding wheat to obtain flour.
- By analyzing how the wheat performs in the experimental mills, this gives a good image of how it will perform in the industry.
- Based on this, it allows select wheats or blend them not only looking at extraction rate but considering the complete behavior during milling.











